




# SEILER

MEDICAL DIVISION  
DENTAL MICROSCOPE



**"YOU CAN'T TREAT  
WHAT YOU CAN'T SEE"**



- 6 step turret magnification system (3 step turret available)
- LED Illumination + 150,000 LUX, 50,000 hour life
- Widefield eyepieces with dioptic locks: 10x, 12.5x 16x, and 20x
- True Coaxial Illumination="Shadow Free"
- Objective Lens: 175mm, 200mm, 250mm, 300mm, and 400mm  
OPTIONAL: Vario Focus Lens - 200mm-350mm
- 0-220° Inclinable Binocular Head (Fixed 45° head available)
- Composite and Red Free Filter standard
- USB and 12V DC power support directly on optic pod (charge Handycams, GoPro, DSLR, HD CCD)

# ALPHA AIR SERIES

FOCUS ON PRECISION.



# DENTAL OPERATING MICROSCOPE

## ACCESSORIES



ASSISTANT  
TEACHING HEAD



0-220°  
BINOCULAR HEAD



WIDE-FIELD  
EYEPIECES



LASER FILTERS



BEAMSPLITTER



EXTENDER  
ROTATION DEVICE



IRIS DIAPHRAGM



GO PRO HERO



DSLR CAMERA  
ADAPTER



HD LIVE VIDEO  
CAMERA



S-VIDEO CAMERA



VARIO FOCUS

## MOUNTING OPTIONS



FLOOR



HIGH WALL



WALL



CEILING



TABLETOP

