

# COMPOUND MICROSCOPES

BY SEILER MEDICAL DIVISION



## MICROLUX IV (MLX817)

### Specifications:

- 5 Infinity Corrected Plan Objectives: 4x, 10x, 20x, 40x, and 100x oil
- Seidentopf Binocular Head, 30° Inclined, 360° Rotation
- Kohler Illumination System: Halogen 12V-20W (LED optional)
- Eyepiece: WF10X/20mm
- Interpupillary Distance: 50-75mm
- Coaxial Coarse and Fine Focusing Control
- Nosepiece: Reversed Quintuple Click-Stop Revolving Mechanism with Multiple Ball Bearings
- Condenser: Abbe Type N.A.=1.25 With Iris Diaphragm
- Double Layer Mechanical Stage with Removable Slide Holder, 190x140mm
- Travel Stage Area 78x55mm with A Right-Hand Stage Handle
- Anti-Fungal Systems
- 110V-240V Universal Power Supply
- Limited Lifetime Warranty

Optional video or digital configurations as well as phase contrast, LED, and MOHS for dermatology procedures.

## WESTLAB III (WLB814)

### Specifications:

- Infinity-Plan Objectives: 4x, 10x, 40x, and 100x oil
- 12V 20W Halogen Illumination
- Seidentopf Binocular Head, 30° Inclined, 360° Rotation
- Eyepiece: WF10X/20mm
- Interpupillary distance: 50-75mm
- Reversed Quadruple click-stop nosepiece
- Coaxial Coarse And Calibrated Fine Focus Control
- Condenser: Abbe Type N.A.=1.25 With Iris Diaphragm
- Mechanical stage with removable slide holder 142 x 140mm; travel area with a right-hand stage handle
- Anti-Fungal System
- 100V – 240V universal power supply
- Optional: Live Video Feature
- Limited Lifetime Warranty

Optional video or digital configurations as well as phase contrast.

## SEILERSCOPE (SXS820)

### Specifications:

- Achromat Objectives: 4x, 10x, 40x, and 100x oil
- Seidentopf Binocular Head, 30° Inclined, 360° Rotation
- Eyepieces: WF10X/18mm
- 6V20W Halogen Illumination
- Coaxial Coarse And Calibrated Fine Focus Control
- Graduated, Right Hand Double Plated Mechanical Stage
- Reversed Four-Hole Nosepiece
- Focusing Tension Adjustment Knob
- Condenser: Abbe Type N.A.=1.25 With Iris Diaphragm
- 100V – 240V universal power supply
- Anti-Fungal System
- Optional Accessories: Live Video, Plan Objectives
- Limited Lifetime Warranty

Optional video or digital configurations as well as phase contrast.

