



SEILER

MEDICAL DIVISION

COMPOUND MICROSCOPES



COMPOUND MICROSCOPES

MICROLUX IV SERIES



Specifications:

- 5 Infinity Corrected Plan Objectives: 4x, 10x, 20x, 40x and 100x Oil (2.5x, 50x and 60x available)
- Seidentopf Binocular Head, 30° Inclined, 360° Rotation
- Kohler Illumination System: Halogen 12V-20W (LED optional)
- Eyepiece: WF10X/20mm
- Interpupillary Distance: 50-75mm
- Coaxial Coarse and Fine Focusing Control
- Nosepiece: Reversed Quintuple Click-Stop Revolving Mechanism with Multiple Ball Bearings
- Condenser: Abbe Type N.A.=1.25 With Iris Diaphragm
- Double Layer Mechanical Stage with Removable Slide Holder, 190x140mm
- Travel Stage Area 78x55mm with Right-Hand Stage Handle
- Antifungal Systems
- 110V-240V Universal Power Supply
- Limited Lifetime Warranty

Available Configurations:

- Camera Attachments - DSLR, HD Video and CCD Video
- Dual Teaching Head
- MOHS for Dermatology Procedures (2.5x, 4x, 10x, 20x, 40x)
- Phase Contrast
- Simple Polarization Kit



WESTLAB III SERIES

Specifications:

- Infinity Plan Objectives: 4x, 10x, 40x and 100x Oil (20x and 50x oil option)
- 12V 20W Halogen Illumination
- Seidentopf Binocular Head, 30° Inclined, 360° Rotation
- Eyepiece: WF10X/20mm
- Interpupillary Distance: 50-75mm
- Reversed Quadruple Click-Stop Nosepiece
- Coaxial Coarse And Calibrated Fine Focus Control
- Condenser: Abbe Type N.A.=1.25 With Iris Diaphragm
- Mechanical Stage with Removable Slide Holder 142 x 140mm; Travel Area with a Right-Hand Stage Handle
- Antifungal System
- 100V – 240V Universal Power Supply
- Optional: Live Video Feature
- Limited Lifetime Warranty

Available Configurations:

- Camera Attachments - DSLR, HD Video and CCD Video
- Dual Teaching Head
- Simple Polarization Kit

SEILERSCOPE SERIES

Specifications:

- Achromat Objectives: 4x, 10x, 40x and 100x Oil
- Seidentopf Binocular Head, 30° Inclined, 360° Rotation
- Eyepieces: WF10X/18mm
- 6V20W Halogen Illumination
- Coaxial Coarse and Calibrated Fine Focus Control
- Graduated, Right-Hand, Double-Plated Mechanical Stage
- Reversed Four-Hole Nosepiece
- Focusing Tension Adjustment Knob
- Condenser: Abbe Type N.A.=1.25 With Iris Diaphragm
- 100V – 240V universal power supply
- Antifungal System
- Optional Accessories: Live Video, Plan Objectives
- Limited Lifetime Warranty

Available Configurations:

- Camera Attachments - HD Video and CCD Video
- Plan objectives available

Seiler Medical Division
3433 Tree Court Industrial Blvd.
St. Louis, MO 63122

1-800-489-2282
1-314-968-2282

www.seilermicro.com



FM 680767



www.seilermicro.com

MED 8054 Rev.0 Sept. 2018