

# NEW CO<sub>2</sub> LASER BARRIER FILTER



The Seiler Laser Barrier Filter (Part # Z-LBFD-2) Device is intended for use with the Seiler Dental or Medical microscopes.

This laser barrier filter was designed to work with the CO<sub>2</sub> Lasers. This small compact filter causes no color distortion and no change in visuals when placed on the microscope, thus providing superior visualization.

Protects against laser wavelengths ranging from:

OD 7@ 9000-11000nm

OD .49@ 520-532nm

Protects against UV light (shade 3) plume ranging from:

180-200nm <0.1% transmission

200-380nm <3.7% transmission



Contact Seiler for more information

1-800-489-2282 • [micro@seilerinst.com](mailto:micro@seilerinst.com) • [www.seilermicro.com](http://www.seilermicro.com)

